USCI

US COMPUTERS, INC.

Established in 1997, a Disadvantaged, Woman-Owned, Small Business









U S Computers, Inc., is an authorized distributor of REHADAPT and offers its entire line of customized products and related items through the ESC contract.

These products are specialized for people with movement impairment.


US Computers, Inc.

PHONE (800)550-8783: FAX (714)528-0530

INFO@USCOMPUTERSING.COM

181 WEST ORANGETHORPE AVE SUITE C PLACENTIA, CA 92870






Unimpeded communication

Rehadapt tries to support dedicated programs to help individuals with special communication needs as often as possible.

In the course of celebrating Saltillo Corporation's jubilee – a manufacturer of communication devices in the AAC field for 20 years now – Rehadapt was honored to present an individual mounting solution for the jubilee winner Alex new communication device – customized and perfectly adapted through our unique Virtual Mounting Solution (VMS) service.

"At school we are excited to use the mount because no matter which classroom he is in the NovaChat can be in the right position for his visual and accessibility needs. It fits close to the chair so he can move around tight classrooms and he can even pull his chair up to the classroom tables."

Jesus Alejandro Reyes Sanchez (Alex) with his teacher Lisa Peterson





Sara Holjević is 8 years old, lives in Karlovac and has cerebral palsy. She can not walk or talk but communicate. Using a speech generating device mounted by Rehadapt Sara greatly improves communication with the environment, can clearly express her feelings, needs and desires, to learn, to answer questions and play games. This fall she is moving into the first grade of school and we believe that our mounting solution is a great help in learning and mastering the program, and the successful establishment of communication in the new environment


Sara Holjević with her parents

WHEELCHAIR MOUNTS

There a is a plethora of wheelchairs being used, with many significant differences. We get your device securely mounted on any of them! A huge range of base clamps ensure solid fixation to any manual or powered wheelchair: clamping onto the frame, on the arm- or footrest, or underneath the seat. Our system is easy to use and quick to adjust even without tools.



Three-tube solutions for most common situations



Standard

Double tubes for more rigidity at higher load



Eye Control Short mounting length ideal for children's wheelchairs



Short

Folding mount for easy exiting and parking



Fold down For extra offset without using additional joints

5


Curved

One- to threetube-solutions for lightweight

devices



Standard



EyeControl

EyeControl mounts include longer tubes to allow for the positioning of your communication aid precisely where you need it. Available in two versions for devices lighter or heavier than 3,2 kg (7 lbs).



strut for more rigidity with higher loads. Perfect solution designed for eye tracking and other devices heavier than 3.2 kg (7 lbs).

TABLE MOUNTS

Our tablemount product range ensures very sturdy mounting for all types of devices. Whether you want to place the device on top of something – even on uneven surfaces – or whether you want to clamp it onto a piece of furniture, Rehadapt has the right solution for you!

Tablestands Self-standing tablemounts that hold devices on

Self-standing tablemounts that hold devices on a tripod, in a kickstand or on an x-frame which automatically brings the mount into an upright position on uneven surfaces.

Tableclamps Tableclamps are securely clamped onto a

Tableclamps are securely clamped onto a given table or armrest. This saves space on the surface. Optionally they can be mounted on walls.



L3D

Light 3D tablemounts with one, two or three tubes each equipped with either joints with Allen key screws (AK) or a QuickShift levers (QS).

Available with either the universal device socket (UDS) and spigot link system (SLS). Combine with any SLS/UDS cradle or device adapter. All equipped with the SuperClamp.



	SLS		UDS	
	Allen screw	QuickShift	Allen screw	QuickShift
1 tube	L3D-TC 1AK	L3D-TC 1QS	L3D-TC 1AK-U	L3D-TC 1QS-U
	14.4001	14.4051	14.4101	14.4151
2 tubes	L3D-TC 2AK	L3D-TC 2QS	L3D-TC 2AK-U	L3D-TC 2QS-U
	14.4002	14.4052	14.4102	14.4152
3 tubes	L3D-TC 3AK	L3D-TC 3QS	L3D-TC 3AK-U	L3D-TC 3QS-U
	14.4003	14.4053	14.4103	14.4153





Anne

Anne Sophie is happy with her Light 3D mount from Rehadapt. She told us it's not just practical but it's the only mount which is also looking cool.



Anne Sophie happy with her new L3D mount from Rehadapt



Wall mounts

Wall mounting plate with base for mounting Monty 3D tubes. Fixation screws and dowels not included.



"M3D Wall mount Base Perpendicular" for mounting tubes perpendicularly to the wall.

11.1051



"M3D Wall mount Base Parallel" mounting tubes parallel to the wall. Tube clampes with a QuickShift lever.

11.1052



"M3D Storage Wall mount" This wallmount unit consists of a wall mounting plate with base and a sturdy M3D tube. The mounting bracket allows Monty 3D wheelchair mounts to be removed from the wheelchair and packed away on the wall when their usage on the wheelchair is not needed.

11.5001



"M3D Wall mount"

Wall mount that includes two tubes (one of them L-shaped), an universal joint and an universal device socket (UDS). Includes slipout safeguards on all tubes and QuickShift levers on all joints. Must be mounted on stable walls. Fixation screws and dowels not included.

11.5000

Trainings

To position a communication device sounds as challanging to you as climbing the Mount Everest? That is what we want to change!


"How to mount an eyetracker to a powered wheel-chair full of special controls and mechanical features?" "What is the best location for a mount to clamp?" "How can I add a second switchmount without wasting space for other wheelchair equipment?"

Together we will find a way!

While we try to simplify the work of our partners with our VMS service and our new service kit, achieving the perfect mount can still be challenging. The best way to become an expert is through hands-on experience. That is why Rehadapt started organizing regular training events - the Rehadapt Mounting Days.

In these hands-on workshops we elaborate together to find the quickest and most expedient mounting solution. With a wide range of wheelchair types you will learn real-world situations as you find in the field. The course will allow you to broaden your knowledge and adding efficiency, economy and quality in providing mounting solutions.

If you are interested and want to know when the next Rehadapt Mounting Days will start just give us a call!



Rehadapt's Hands-on-Workshops



Adam, his family, teachers and therapists are so excited about his new communication device from Tobii Dynavox and the mounting arm from Rehadapt Engineering that will mean he can always have it with him wherever goes. Thank you to the teams at both of these remarkable companies that made this possible.

Inclusive Solutions, South Africa



Adam with his teachers and therapists



One of the users of our mounts, 10 year old Julia Nickel, became a celebrity after a video showcasing her and other Assistive Technology users gathered more than 1.5 million (!) views on Youtube. When she met us at Rehadapt's booth at Rehacare, we serviced her mount. In the following panel discussion she gave insights into Alternative and Augmentative Communication to a general audience.

We are glad you stopped by, Julia!

Julia Nickel with Rehadapt Engineer Marco Miftari

Boston Children's Hospital

"We've LOVED having the flexibility of the Rehadapt rolling floor mount at bedside for our inpatient AAC service. The ability to use systems such as the Tobii i12 device pictured here (in a black case) in multiple positions at bedside with a very fast, easy and reliable position adjustment that conforms to the child's varied positions in bed has been great!"

John Costello, MA, CCC-SLP Director of the Augmentative Communication Program at Boston Children's Hospital with Uli Ehlert, CEO Rehadapt Engineering



What is a good **mount** without a great **SETVICE**



Our unique service: we analyze your specific needs to provide you with the best mounting solution - **for free!**

The way to your individual Virtual Mounting Solution:

1. Ask for it!

Send us pictures of your wheelchair and basic information about the chair and your device

2. We engineer

Our experienced engineers will design a custom mounting solution within just three business days

- 3. After 3 days we get back to you
 We send you a proposal with the visualized solution photoshopped into your picture
- 4. And we deliver after 5 days!
 We ship your pre-assembled individual solution within five business days after order



✓ Fast

only 3 business days and you can see a visual proposal in detailed superimposed images

✓ Economical

no meetings, no travelling efforts, no wait times required ✓ Simple

we provide illustrated assembly instructions and an easy-to-handle order ID

Start right now!









Your personal mount customized by our engineers to meet your individual needs



What we offer to you —
we analyze your specific needs
to provide you the best mounting solution — for free!



Fast

only 3 business days and you can see a visual proposal in detailed superimposed images



no meetings,
no travelling efforts,
no wait times
required



Simple
we provide
illustrated assembly
instructions and an
easy-to-handle order ID





The way to your personal solution



You ask for it

Send us pictures of your wheelchair and basic information about the chair and your device

We engineer...

Our experienced engineers will design a custom mounting solution within just three business days

...get back to you...

We send you a proposal with the visualized solution photoshopped into your picture

...and deliver

We ship your pre-assembled individual solution within five business days after order

You win!



A unique VMS ID allows for easy order handling and after sales service. In case of changing the wheelchair or device, full traceability saves time and money

Trust the global leader of wheelchair mounts.

Engineering the best mounting solution for our customers is our passion. You will benefit from our devotion to quality during long years of usage. All our products are manufactured in Germany, using high tenacity Aluminium alloys as in leading-edge wheelchairs and aircraft design.

Extensive knowledge, high-tech materials and ultra-precise manufacturing are the components to build what sounds so simple:

The perfect mount for you. Lightweight, rugged and personal.


More 7 benefits.





US COMPUTERS, INC.

EXPERIENCE THE PROFICIENCY OF OUR SERVICE TODAY



Established in 1997. A Certified Small Business with State of California

Our Contracts ensures that government agencies can purchase our products and services with confidence and assurance that we are a reliable company offering our goods and services at competitive prices approved.

Woman-owned, disadvantaged, small business

US Computers, Inc.

(800)550-8783